

# NE-1 Neon Calibration Light Source Installation and Operation Instructions

## Description

The NE-1 Neon Calibration Source is a wavelength calibration source for UV-VIS-Shortwave NIR spectrophotometric systems. The NE-1 produces Neon lines from 540-754 nm, and is an ideal lamp to use when performing fast, accurate spectrometer wavelength calibrations.

The NE-1 features an SMA 905 Connector for interfacing with our optical fibers. It operates with a 12 VDC power supply (included with the unit) or 9V battery (not included).

The following sections detail the features of the NE-1 Neon Calibration Light Sources.



**NE-1 Neon Calibration Light Source**

---

### Note

The NE-1 is NOT designed to operate as an excitation source in your experiments. Spectral lines above 922 nm are not easily detectable by Ocean Optics spectrometers.

---

# Parts Included

The NE-1 package includes the following items:

- NE-1 Neon Calibration Light Source
- 12 VDC power supply

---

## WARNING

**The beam emerging from the light source contains UV radiation that can cause serious eye injury upon direct contact with the eye. Never look directly into the light source.**

**The SMA 905 Connector may get extremely hot during operation. After lamp use, allow sufficient time to cool before handling.**

**Dangerous voltages are present, and there are no user-serviceable parts inside. Never open the NE-1.**

---

## Additional Accessories

The following are additional accessories available from Ocean Optics that you may need, depending on your system set-up:

- Spectrometer
- SMA-terminated optical fiber
- Ocean Optics software

## Connecting the NE-1

Follow the steps below to set up your NE-1 for use.

### Procedure

1. Plug the 12 VDC power supply into a power outlet, then connect the barrel connector of the power supply to the power input on the rear of the NE-1.

or

Alternately, you can use a 9-volt battery (not included) to power the NE-1. Open the battery hatch of the NE-1 and install the 9-volt battery, and then proceed to Step 2.

2. Connect a fiber to the SMA 905 Connector on the NE-1. If your spectrometer does not have an entrance slit, use a 50  $\mu\text{m}$  diameter (or smaller) fiber. Larger fibers and slits result in reduced optical resolution.

### Note

If the spectrometer does not have a slit and your experiment requires you to use fibers of varying diameters, you will need to perform a wavelength calibration after changing fibers. You should perform a wavelength calibration each time you unscrew the fiber from the spectrometer and change fiber size.

3. Move the On/Off switch on the NE-1 (next to the SMA 905 Connector) to the On position. The red LED will illuminate to indicate that the NE-1 is powered on.

You have now configured the NE-1 for use.

## Calibrating With the Light Source

The information in this section explains how to calibrate your spectrometer's wavelength using the NE-1 light source.

### About the Wavelength Calibration

You are going to be solving the following equation, which shows that the relationship between pixel number and wavelength is a third-order polynomial.

$$\lambda_p = I + C_1p + C_2p^2 + C_3p^3$$

Where  $\lambda_p$  = the wavelength of pixel  $p$ ,  
 $I$  = the wavelength of pixel 0,  
 $C_1$  = the first coefficient (nm/pixel),  
 $C_2$  = the second coefficient (nm/pixel<sup>2</sup>)  
 $C_3$  = the third coefficient (nm/pixel<sup>3</sup>)

You will be calculating the value for  $I$  and the three  $C$ s.

## Calibration Requirements for Ocean Optics Spectrometers

To re-calibrate the wavelength of your Ocean Optics spectrometer using the NE-1, you will need the following items:

- NE-1 Neon Calibration Light Source
- Ocean Optics spectrometer and its manual
- An optical fiber (for spectrometers without a built-in slit, a 50- $\mu$ m fiber works best)
- A spreadsheet program (Excel or Quattro Pro, for example) or a calculator that performs third-order linear regressions

**Note**

If you are using Microsoft Excel, choose **Tools | Add-Ins** and check **AnalysisToolPak** and **AnalysisToolPak-VBA**.


## Calibration Procedure

Please see the appropriate manual for your model Ocean Optics spectrometer for instructions on how to calibrate it. Ocean Optics manuals are located at <http://oceanoptics.com/support/technical-documents/>

## NE-1 Specifications

The following sections detail the specifications of the NE-1.

### Spectral Output



## Strong NE Emission Lines by Wavelength (nm)

There are more emission lines shown below than printed on the label on the NE-1 housing. The label is intended as a quick, convenient reference and does not list every NE emission line that exists.

Below is a list of wavelengths for Neon (as well as Krypton and Xenon). Wavelengths listed in blue are the ones listed on the calibration lamp's label.

| Wavelength | Lamp | Wavelength | Lamp | Wavelength | Lamp | Wavelength | Lamp |
|------------|------|------------|------|------------|------|------------|------|
| 341.790    | Ne   | 503.135    | Ne   | 733.930    | Xe   | 904.545    | Xe   |
| 342.391    | Ne   | 503.775    | Ne   | 738.600    | Xe   | 916.265    | Xe   |
| 344.770    | Ne   | 508.038    | Ne   | 739.379    | Xe   | 979.970    | Xe   |
| 345.076    | Ne   | 511.367    | Ne   | 740.040    | Xe   | 992.319    | Xe   |
| 345.419    | Ne   | 511.650    | Ne   | 743.890    | Ne   | 1083.837   | Xe   |
| 346.052    | Ne   | 540.056    | Ne   | 747.244    | Ne   | 1117.752   | Ne   |
| 346.658    | Ne   | 556.222    | Kr   | 748.887    | Ne   | 1152.275   | Ne   |
| 347.257    | Ne   | 557.029    | Kr   | 753.577    | Ne   | 1181.938   | Kr   |
| 349.806    | Ne   | 576.441    | Ne   | 754.404    | Ne   | 1220.353   | Kr   |
| 350.121    | Ne   | 582.015    | Ne   | 755.979    | Xe   | 1317.741   | Kr   |
| 351.519    | Ne   | 585.249    | Ne   | 758.468    | Xe   | 1363.422   | Kr   |
| 352.047    | Ne   | 587.096    | Kr   | 758.741    | Kr   | 1365.706   | Xe   |
| 359.353    | Ne   | 588.189    | Ne   | 760.155    | Kr   | 1414.244   | Xe   |
| 360.017    | Ne   | 594.483    | Ne   | 764.391    | Xe   | 1442.679   | Kr   |
| 363.366    | Ne   | 597.553    | Ne   | 768.525    | Kr   | 1473.281   | Xe   |
| 368.573    | Ne   | 602.999    | Ne   | 769.454    | Kr   | 1473.444   | Kr   |
| 370.122    | Ne   | 607.433    | Ne   | 780.265    | Xe   | 1537.204   | Kr   |
| 431.958    | Kr   | 609.616    | Ne   | 785.482    | Kr   | 1541.839   | Xe   |
| 436.264    | Kr   | 612.884    | Ne   | 788.132    | Xe   | 1605.328   | Xe   |
| 437.612    | Kr   | 614.306    | Ne   | 791.343    | Kr   | 1672.815   | Xe   |
| 439.997    | Kr   | 616.359    | Ne   | 796.734    | Xe   | 1689.676   | Kr   |
| 445.392    | Kr   | 621.728    | Ne   | 805.726    | Xe   | 1800.223   | Kr   |
| 446.369    | Kr   | 626.649    | Ne   | 805.950    | Kr   | 1816.733   | Kr   |
| 450.235    | Kr   | 630.479    | Ne   | 806.134    | Xe   | 2190.851   | Kr   |
| 452.186    | Xe   | 633.442    | Ne   | 810.436    | Kr   | 1262.339*  | Xe   |
| 462.420    | Xe   | 638.299    | Ne   | 819.006    | Kr   | 1520.310*  | Kr   |
| 466.849    | Xe   | 640.225    | Ne   | 823.163    | Xe   | 1620.872*  | Kr   |
| 469.097    | Xe   | 650.653    | Ne   | 826.324    | Kr   | 1647.29*   | Xe   |
| 469.804    | Xe   | 653.288    | Ne   | 826.652    | Xe   | 1656.023*  | Xe   |
| 473.415    | Xe   | 659.895    | Ne   | 829.811    | Kr   | 1755.350*  | Kr   |
| 479.262    | Xe   | 667.828    | Ne   | 837.761    | Ne   | 1763.882*  | Xe   |
| 480.702    | Xe   | 671.704    | Ne   | 849.536    | Ne   | 1785.738*  | Kr   |
| 482.971    | Xe   | 692.947    | Ne   | 877.675    | Kr   | 1790.45*   | Xe   |
| 484.329    | Xe   | 703.241    | Ne   | 878.375    | Ne   | 1809.09*   | Xe   |
| 491.651    | Xe   | 717.394    | Ne   | 881.941    | Xe   | 1832.53*   | Xe   |
| 492.315    | Xe   | 724.512    | Ne   | 892.869    | Kr   | 1959.94*   | Xe   |
|            |      |            |      |            |      | 1984.638*  | Xe   |

**Specifications Table**

|                          |  |
|--------------------------|--|
| Output                   | Low-pressure gas discharge lines of Neon                       |
| Dimensions (in mm):      | 125.7 x 70 x 25.8  |
| Power consumption:       | 250 mA at 12 VDC   |
| Power requirements:      | 12 VDC wall transformer (included) or 9 VDC battery (optional) |
| Bulb life:               | Approx. 3500 hours (at 20 mA)                                  |
| Internal voltage:        | 600 volts at 30 kHz  |
| Aperture:                | 3 mm   |
| Amplitude stabilization: | ~ 1 minute   |
| Connector:               | SMA 905  |